



Dear Patient,

the monitoring of your personal values for blood coagulation by an experienced physician is for your health of high importance. The transmission of those data by mobile communication networks makes you sure that the latest data always are evaluated by a medical specialist. This unique technology is very comfortable for you and can be used all over the world, whether you are at home or travelling.

This is how you transmit the health data into your personal data base:

1. Check the values for blood coagulation as you are usually doing with the **CoaguChek - XS®** meter.
2. To transfer the health data, please put the **CoaguChek - XS®** with the front side in direction to the Life Data Agent.



Switch on the **CoaguChek - XS®** by pressing the on/off- button.

3. Switch on the Life Data Agent by keep pressed the **ON- Button** until the device starts (LED are blinking 3 times).

Notice: **CoaguChek - XS®** is a registered trademark of Roche Diagnostics



4. About 10 seconds after switched on the **first LED** (left) starts blinking thus indicating that the **Life Data Agent** is searching the **mobile communication network**. The LED will stay lighted when this connection is established.



5. Blinking of the **second LED** (middle) shows that the **connection to the data base** is in progress. The LED stay lighted when the connection to the data base is established.



6. The next step is the **transmission of your health data**. The Life Data Agent transmits all data automatically. This process is indicated by the **third LED** (right).



7. As soon as all data are transmitted successfully **all three LED** are switched on permanently. A few seconds later the Life Data Agent switches off automatically.



8. If you want switch off the Life Data Agent manually (e.g. if you are in a hurry), simply press the **Off button**.

**Accu charge:** The Life Data Agent was given to you with a fully charged Li-Polymer-accu. After about two weeks of use you should charge the accu. A preventive charging has no negative impact and is useful as well.

For charging please use the power supply (wall adaptor) that is delivered with the Life Data Agent. The LED at the right side of the device is indicating the charging process: This LED is **permanently on** when charging. It starts **blinking** when the accu is full and the charge process performs the trickle-charge mode in order to keep the accu charged.



## If something does not work.....

### **1. The device is starting as usual, but the first LED (left) does not work permanently**

#### Reason for the fault:

The Life Data Agent has no connection to the cellular network, because the performance of the antenna is bad.

#### How to manage:

Try it again later and move to a place where a stable mobile communication is possible.

### **2. The device is starting, but the LED in the middle stays blinking all the time.**

#### Reason for the fault:

The device has contact to the mobile communication network, but is unable to connect to the data base. Maybe the data base is in a service routine.

#### How to manage:

Please try later to transmit your data.

### **3. All LED are glowing weakly**

#### Reason for the fault:

The accuracy of the Life Data Agent is too weak and the device is unable to send data.

#### How to manage:

Please charge the accu and then try again to transmit your health data.